

# CURESp<sup>ON</sup>SE

**Curesponse is  
on a Mission to  
Improve Cancer  
Treatment.**



**cResponse™ combines rapid genomic  
sequencing with a functional test,  
enabling prioritization of different therapeutic  
options for patients with metastatic cancer.**

# The cResponse™ solution


We offer an innovative functional test that in combination with rapid genomic sequencing, could assist your oncologist in prioritizing treatments and selecting the most effective treatment for you.

## Background

It is well established that all tumors are different, and each patient will respond differently to anti-cancer drugs. Genomic sequencing provides a map of the genetic profile of your tumor and the potential therapeutic options. The cResponse™ test actually evaluates them on your tumor and assesses their relative effect.

The tumor microenvironment (TME), plays a central role in drug response and in many cases is a major factor in developing resistance of cancer to drugs. The cResponse™ platform preserves the tissue's 3D structure and maintains the TME, enabling assessment of your tumor's response to drugs in laboratory conditions.

## The process





A fresh biopsy arrives at a regional cResponse™ service center. Part of the sample undergoes rapid genomic sequencing to identify tumor-specific genomic alterations that may be targeted by specific anti-cancer drugs. Added to those selected by the physician, they create a panel of therapies to be tested in the system. The biopsy is cut into thick sections which are placed in a proprietary assay culture matrix and tests the cancer's sensitivity to various drugs or novel drug combinations. cResponse's proprietary evaluation system and response algorithm then prioritizes the options according to demonstrated drug efficacy, providing results, usually within two weeks.



# The report


Upon conclusion of the test and analysis, we will provide a comprehensive report, summarizing the genomic information found in our rapid 57-gene panel, the different drugs we tested and the effect of each tested drug on the tumor sample as assessed in our platform.



The cResponse™ test is using a proprietary algorithm, where each drug is assigned a score of 0-100, with 100 reflecting a strong response to therapy in our assay and 0 reflects no response. The highest scoring drugs represent the most potent treatment option as tested in our proprietary assay and may indicate higher chances of clinical effect.

## Drug response chart


### cResponse™ Score





Images of the cancer tissue are also included in the report, demonstrating tissue morphology, viability and replication.

## Images of the cancer tissue under different treatments



## A list of the genes tested in our panel

### 286 amplicons covering hotspots of 57 actionable genes

ABL1, AKT1, ALK, APC, ATM, BRAF, CDH1, CDKN2A, CSF1R, CTNNB1, DDR2, DNMT3A, EGFR, ERBB2, ERBB4, EZH2, FBXW7, FGFR1, FGFR2, FGFR3, FLT3, FOXL2, GNA11, GNAQ, GNAS, HNF1A, HRAS, IDH1, IDH2, JAK2, JAK3, KDR, KIT, KRAS, MAP2K1, MET, MLH1, MPL, MSH6, NOTCH1, NMP1, NRAS, PDGFRA, PIK3CA, PTEN, PTPN11, RB1, RET, SMAD4, SMARCB1, SMO, SRC, STK11, TP53 (Whole CDS coverage), TSC1, TSC2, VHL.

[support@curesponse.tech](mailto:support@curesponse.tech)

### Curesponse:

2 Oppenheimer street,  
Haogen Building B, 14th floor,  
Rehovot, Israel  
Tel: 972-74-764610

### Curesponse UK Ltd:

Translation & Innovation Hub  
Lab 1 19 L, Imperial College  
White City Campus, 84 Wood Lane  
London W12 0BZ